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
<http://www.dentwise.eu/about/viewer>

### 2. General Functions

Start the application by double clicking on the icon on the desktop or via the start menu > All Applications > ATOS Viewer v#.#.# > ATOS Viewer

On the left side of the application you will find the *explorer* showing the different meshes. The viewer sent to you contains **3 meshes: tooth preparation mesh (\*\_TO), gums mesh (\*\_TV) and a mesh containing the designed framework (\*\_Viewer).**

Use the *toggles* to show or hide a mesh:

The right half of the application shows a *3D view* showing the visible meshes.


Use the **left mouse button** to **rotate the 3D view**.

**Scroll** the mouse wheel to **zoom the 3D view**.

Click the **mouse wheel** or **middle mouse button** to **pan the 3D view**.

The lower left corner of the application window shows the *Info* tab and the *3D Overview* tab. The 3D overview tab allows you to [make section views](#).

### 3. Measure Distances

The viewer enables you to measure distances in various ways via the menu *Distances*.


For simple **measurements between two points**, use *Distances > Point-Point Distance*. Select a point by pressing the **CONTROL** button and **left mouse click** on the mesh. Do the same for Point 2. You can change Point 1 and 2 by selecting a point again until you have created the distance. **Hide meshes** you don't need and **rotate the 3D view** to be able to select the point that you want. **Click Create** to keep the measurement.

Measurements appear in the explorer under *Dimensions > Distances*. The toggle enables you to show or hide measurements individually.

To **remove elements (dimensions, meshes...)** from the explorer, right click on the element and choose *Close* in the submenu.

## 4. Make Section Views

Enable the **3D Overview** by clicking 3D Overview tab in the left bottom corner of the application window. **Right click** in the 3D overview shows this submenu:


By selecting the X or Y plane, you can cut the meshes in the frontal or distal area to check the available space.<sup>2</sup>

Right click in the 3D overview *Visibility > Visibility from Explorer* allows you to **show or hide meshes in the 3D overview**.

<sup>2</sup> To make your own section view: click **User-Defined Clipping Plane** in the submenu. Click in the 3D overview to choose a start point and **draw a line** by clicking again as end point of the line. A section clipping plane will be created through the line, perpendicular to the 3D overview screen.

## 5. Changing Mesh Visibility Options

Right click on a mesh in the explorer and choose *Edit Properties* to change the appearance of mesh:


Mesh rendering:

Points rendering:

Triangles rendering:

Click the *Rendering* tab: here you can display the mesh as points, mesh or triangles, or a combination of these. Here you can also change colors, point size...

## 6. Save Your Work

You can **save the session** and review your changes at a later time or send the file to someone else: menu *File > Save Session* or *File > Save Session As*.

## 7. Close Session

**Close a complete session** via the menu *File > Close All*.

Contact us if you need more help or information.